


ApolloAnkle™

Fracture Plating System



Seeing is **Believing**



Tapered proximal end easily and comfortably fits under soft tissue upon insertion

Special syndesmotic holes are designed for the 4.3mm screws and standard suture buttons. They are anatomically angled 35°, posterior to anterior

Screw heads sit below the plate surface when fully engaged

The plates are built with a patented hybrid construction of 3D printed titanium and injection molded PEEK. The combination of these manufacturing techniques, along with a blend of proprietary surface treatments, creates a lightweight, bendable plate with strength equal to traditional machined titanium plates with Type II anodization.

Ø2.9mm and Ø3.7mm locking and non-locking screws can be inserted up to 15° off-axis

Plate holes allow for Ø2.9mm and Ø3.7mm locking and non-locking screws

Proprietary Hybrid Titanium/PEEK construction is ortholucet and malleable

*PEEK pellets are melted and injected into a custom mold. Any left over sprues or scrap can be reground into the raw material and used again. Components produced using additive manufacturing require minimal to no additional machining and create no titanium waste.

Hybrid Materials and Manufacturing

The plates are not only visually green, but also environmentally green. Produced using a zero waste process, these plates are made through a proprietary additive manufacturing (3D printing) and injection molding process.

Additive manufacturing offers the capability to mirror complex anatomy and create very smooth contours that minimize soft tissue irritation. This technology creates the shell structure of the plate which offers the thinnest section of titanium on the market, giving Apollo its radiolucent properties.*



Low-profile, contoured plate offers strength while reducing soft tissue irritation

Distal fibula screw cluster allows for multiple points of fixation

ZERO
WASTE
MANUFACTURING


Screw Technology

The unique combination of additive manufactured titanium and injection molded PEEK leads to the patented screw locking technology of **PEEKLOC™**.

The PEEK construct allows for a quick thread engagement, and smooth tactile feel, while the proprietary hidden titanium structure enforces the solid locking strength. PEEKLOC™ technology creates the elastic locking friction and greatly reduces the risk of cold-welding during screw insertion.

Hybrid Cortical / Cancellous Screw System

The optimized thread design and slightly larger OD provide greater pull out strength in both bone types. So, only one screw type needed, simplifying the surgery.



PEEKLOC™ and Self-countersinking Screws

Below surface insertion of screw

T15 Hexalobe socket

2mm Driver depth

Patented PEEKLOC™ system creates elastic friction when locking screw in place

Self-countersinking features add more depth below the plate for a flush mounted screw with deeper driver engagement

All Apollo locking and non-locking screws can be placed up to 15° off-axis.





PEEKLOC™ provides a quick thread engagement and smooth tactile feel when locking screw in place

PEEKLOC™ creates a strong interface between the screw and plate, with better cantilever strength and accepting higher torque strength

PEEKLOC™ minimizes titanium to titanium interface therefore reducing cold-welding

Partially-cannulated for ortholucency

Wider outside thread diameter gives a hybrid cortical and cancellous screw design




Ortho

This thin cross-sectional structure of titanium gives the plates "ortholucency", a unique advantage over traditional metal implants where the radiolucent properties drastically improve the visualization of bones and joint spaces. Intraoperatively, surgeons benefit from improved fracture and joint reduction imaging, while postoperatively they can better assess if bone is healing properly, potentially leading to earlier weight bearing decisions.



Apollo Ankle Fracture Plating System

Recent Technology



1. Posterior Tibial Plate

- Left/Right anatomical designs.
- Distal holes angled away from tibiotalar joint.

2. Fibular Plate

- Syndesmotomic holes positioned to aim 35°, posterior to anterior.
- Syndesmotomic holes designed to accept typical suture buttons.
- Multiple points of fixation in the distal cluster.

3. One-third Tubular Plate

- Versatile plate with hole choices from 4 to 12 holes.

4. Medial Malleolar Plate

- Extra thin distal portion to minimize soft tissue irritation.


5. Hook Plate

- Versatile design for both lateral and medial malleoli.
- Compression slot with 2mm compressive action.

6. Syndesmotomic Plate

- 2 and 4 hole designs.
- Syndesmotomic hole designed for typical suture button.

Thin Ti AM shell construction



Contoured to the bone and smooth radius on the edges for less soft tissue irritation



Sterile, Disposable Instrument Kits

Optimize your work flow within the surgical procedure.

Single-use Convenience Kit, including two starter screws



T15 Driver	2
Drill Ø3.7mm / Ø4.3mm – Core	1
Drill Ø2.9mm – Core	1
Drill 3.7 – Lag	1
Olive Wire Assembly	3
Ratcheting Handle, Single Use	1
Drill Guide, Polyaxial/Straight	1
Countersink/Depth Gauge,	1
K-Wires, Ø1.6mm x 150mm	2
Plate Benders	2
Ø2.9mm x 12mm Non-Locking Screw	1
Ø3.7mm x 14mm Non-Locking Screw	1

Lag Drill Guide Kits, 2.9, 3.7 and 4.3



Single-use Hook Plate Kit



Single-use Radiopaque Trials




Optional Reduction Instruments*



*Reusable instruments and sterilization tray.



Portfolio Information

	Plate	Length	Shaft Width	Head Width	Hole Count	Orientation
	Lateral Fibula	67 89 111	9.6	18	9 11 13	Left / Right
	One-Third Tubular	51 62 73 84 95 117 139	9.6	N/A	4 5 6 7 8 10 12	Universal
	Hook	48 59 70 81	10	N/A	3 4 5 6	Universal
	Medial Malleolar	60 71 83	9.6	21	6 7 8	Universal
	Posterior Tibia	48 59	10.6	23	5 6	Left/Right
	Syndesmosis	29 51	10.8	N/A	2 4	Universal
	Material	Ti6Al4V / PEEK				



Ø Size range	2.9mm	3.7mm	4.3mm
Type	Locking and Non-Locking		Non-Locking
Length	8-40mm	10-60mm	25-70mm
Material	Ti6Al4V		
Color	Fuchsia	Green	Blue

Screw and Instrument Caddies





Scan for more
product information

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only. Implants and Instrument kits for treatment of temporary and permanent disabled individuals (patients).

GLW, Inc, GLW Medical and Apollo are trademarks of GLW, Inc.. Innov8ortho is the exclusive distributor of the Apollo Ankle Fracture Plating Implants.

www.glwmed.com

© 2025 GLW, Inc. All rights reserved.
Patents: www.glwmed.com/patents/



Legal Manufacturer:

GLW, Inc.
930 Sylvan Ave, Suite 125
Englewood Cliffs, NJ 07632

innov8ortho

Distributed by:

Innov8ortho, LLC
930 Sylvan Ave, Suite 125
Englewood Cliffs, NJ 07632
custsvc@innov8ortho.com

Ref: APO-BRO-07-25-EN