

Medial Malleolar Plate

REF	Hole	Orientation
F3-0003-003S	6	Universal
F3-0003-004S	7	Universal
F3-0003-005S	8	Universal

Posterior Tibial Plates

REF	Hole	Orientation
F3-1002-003S	5	Left
F3-2002-003S	5	Right
F3-1002-004S	6	Left
F3-2002-004S	6	Right

Distal Lateral Fibula Plates

REF	Hole	Orientation
F3-1001-004S	9	Left
F3-2001-004S	9	Right
F3-1001-006S	11	Left
F3-2001-006S	11	Right
F3-1001-008S	13	Left
F3-2001-008S	13	Right

One Third Tubular Plates

REF	Hole	Orientation
F3-0004-004S	4	Universal
F3-0004-005S	5	Universal
F3-0004-006S	6	Universal
F3-0004-007S	7	Universal
F3-0004-008S	8	Universal
F3-0004-010S	10	Universal
F3-0004-012S	12	Universal

Syndesmosis Plates

REF	Hole	Orientation
F3-0005-002S	2	Universal
F3-0005-004S	4	Universal

Hook Plates

REF	Hole	Orientation
F3-0006-003S	3	Universal
F3-0006-004S	4	Universal
F3-0006-005S	5	Universal
F3-0006-006S	6	Universal

F4-3006-000S Single Use Trials Kit

Ø3.7mm Locking and Non-locking Screws



Locking Ø3.7mm	Non-locking Ø3.7mm	Length
F3-1037-010S	F3-0037-010S	10mm
F3-1037-012S	F3-0037-012S	12mm
F3-1037-014S	F3-0037-014S	14mm
F3-1037-016S	F3-0037-016S	16mm
F3-1037-018S	F3-0037-018S	18mm
F3-1037-020S	F3-0037-020S	20mm
F3-1037-022S	F3-0037-022S	22mm
F3-1037-024S	F3-0037-024S	24mm
F3-1037-026S	F3-0037-026S	26mm
F3-1037-028S	F3-0037-028S	28mm
F3-1037-030S	F3-0037-030S	30mm
F3-1037-032S	F3-0037-032S	32mm
F3-1037-034S	F3-0037-034S	34mm
F3-1037-036S	F3-0037-036S	36mm
F3-1037-038S	F3-0037-038S	38mm
F3-1037-040S	F3-0037-040S	40mm
F3-1037-042S	F3-0037-042S	42.5mm
F3-1037-045S	F3-0037-045S	45mm
F3-1037-047S	F3-0037-047S	47.5mm
F3-1037-050S	F3-0037-050S	50mm
F3-1037-052S	F3-0037-052S	52.5mm
F3-1037-055S	F3-0037-055S	55mm
F3-1037-057S	F3-0037-057S	57.5mm
F3-1037-060S	F3-0037-060S	60mm

Ø2.9mm Locking and Non-locking Screws



Locking Ø2.9mm	Non-locking Ø2.9mm	Length
F3-1029-008S	F3-0029-008S	8mm
F3-1029-010S	F3-0029-010S	10mm
F3-1029-012S	F3-0029-012S	12mm
F3-1029-014S	F3-0029-014S	14mm
F3-1029-016S	F3-0029-016S	16mm
F3-1029-018S	F3-0029-018S	18mm
F3-1029-020S	F3-0029-020S	20mm
F3-1029-022S	F3-0029-022S	22mm
F3-1029-024S	F3-0029-024S	24mm
F3-1029-026S	F3-0029-026S	26mm
F3-1029-028S	F3-0029-028S	28mm
F3-1029-030S	F3-0029-030S	30mm
F3-1029-032S	F3-0029-032S	32.5mm
F3-1029-035S	F3-0029-035S	35mm
F3-1029-037S	F3-0029-037S	37.5mm
F3-1029-040S	F3-0029-040S	40mm

Ø4.3mm Non-locking Screws



Non-locking Ø4.3mm	Length
F3-0043-025S	25mm
F3-0043-027S	27.5mm
F3-0043-030S	30mm
F3-0043-032S	32.5mm
F3-0043-035S	35mm
F3-0043-037S	37.5mm
F3-0043-040S	40mm
F3-0043-042S	42.5mm
F3-0043-045S	45mm
F3-0043-047S	47.5mm
F3-0043-050S	50mm
F3-0043-052S	52.5mm
F3-0043-055S	55mm
F3-0043-057S	57.5mm
F3-0043-060S	60mm
F3-0043-062S	62.5mm
F3-0043-065S	65mm
F3-0043-067S	67.5mm
F3-0043-070S	70mm

Washer

F3-0075-000S	7.5mm Fits all screws
--------------	-----------------------



F4-3002-000S
Single Use Lag Kit,
2.9mm screw



F4-3003-000S
Single Use Lag Kit,
3.7mm screw




F4-3004-000S
Single Use Lag Kit,
4.3mm screw



F4-3005-000S
Single Use Hook Plate
Instrument Kit


Optional Reduction Instruments



	REF	Description		REF	Description
1	F5-9002-000	Tray	8	F5-9009-000	Dental Pick
2	F5-9005-000	Syndesmotic Clamp	9	F5-9024-000	Tissue forcep
3	F5-9022-000	Parallel Pin Distractor	10	F5-9021-000	Hohmann Retractor–15mm
4	F5-9006-000	Verbrugge Clamp	11	F5-9007-000	Hohmann Retractor–8mm
5	F5-9004-000	Lobster Clamp	12	F5-9008-000	Periosteal Elevator
6	F5-9020-000	Reduction Forceps – 8"	13	F5-9010-000	Curette, size 00
7	F5-9003-000	Reduction Forceps – 5"	14	F5-9023-000	Distraction Pins, Tube



Push & Lift



F4-3000-000S



Ø3.7mm Sterile Quick Caddy

Ø2.9mm Sterile Quick Caddy



F4-3014-000S



F4-3013-000S

QTY	Description	QTY	Description
2	Ø3.7mm x 12mm Non-Locking	1	Ø2.9mm x 10mm Non-Locking
2	Ø3.7mm x 14mm Non-Locking	2	Ø2.9mm x 12mm Non-Locking
2	Ø3.7mm x 12mm Locking	1	Ø2.9mm x 14mm Non-Locking
2	Ø3.7mm x 14mm Locking	1	Ø2.9mm x 10mm Locking
		2	Ø2.9mm x 12mm Locking
		1	Ø2.9mm x 14mm Locking



REF	Description
F4-3000-000S	Apollo Convenience Kit, with 2 starter screws

Apollo Convenience Kit Contents

	T15 Self-retaining Driver	2
	Drill Ø3.7mm/Ø4.3mm Screw – Core	1
	Drill Ø2.9mm Screw – Core	1
	Drill Ø3.7mm Screw – Lag	1
	Olive Wire Assembly	3
	Ratcheting Handle, Single use	1
	Drill Guide, Polyaxial/Straight	1
	Countersink/Depth Gauge	1
	K-Wires, Ø1.6mm x 150mm	2
	Plate Bender	2

Starter Screws

	Ø2.9mm x 10mm Locking screw	1
	Ø2.9mm x 12mm Locking screw	1

PEEKLOC™

Variable Angle Locking

The PEEK reinforced underside of the plate creates elastic friction for locking the screw and mount below the plate surface where even angled screws can remain flush to the plate.



- No cold welding between plate and screws.
- Screws can pivot freely by $\pm 15^\circ$ in all directions for optimal positioning.
- Self-countersinking adds more depth below the plate.