

# FusionFrame™






Ring Lock System

## Sales Sheet










# FusionFrame™ at a Glance










## External Supports

1	<b>Full Rings</b>	140mm – 205mm	
2	<b>Half Rings</b>	140mm – 205mm	
-	<b>5/8 Rings*</b>	140mm – 205mm	
3	<b>Foot Plates</b>	140mm – 205mm	
4	<b>Walkie Pucks</b>	Tall and short sizes with removable spacer	




## Connectors

5	<b>Threaded Rods</b>	60mm – 400mm	
6	<b>Quick Connect Struts</b>	Three Lengths	
7	<b>Oblique Supports</b>	One Size	
-	<b>Universal Joint*</b>	One Size	
-	<b>Threaded Sockets*</b>	40mm & 60mm	
-	<b>Plates*</b>	20mm & 30mm	
-	<b>90° Hinge*</b>	One Size	





## Fasteners

8	<b>Nuts</b>	Hexagonal Head with 10mm Span	
-	<b>Counting Nuts*</b>	One Size	
9	<b>Bolts</b>	12mm, 16mm & 20mm	
10	<b>Wire Fixation Bolt</b>	Combination Slotted / Cannulated	
-	<b>Half Pin Fixation Bolt*</b>	Accommodates All Half Pin Sizes	
-	<b>Half Pin Fixation Cubes*</b>	1 – 4 Holes	
11	<b>Male / Female Posts and Hinges</b>	1 – 4 Holes	
-	<b>Washers*</b>	Three Thicknesses, Including Slotted Washer	
-	<b>Conical Washer Couples*</b>	One Size	

## Implants

-	<b>Smooth Wire*</b>	1.8mm Diameter	
-	<b>Wire with Stopper (Olive Wire)*</b>	1.8mm Diameter	
-	<b>Half Pins*</b>	5mm & 6mm Thread Diameters	

## Instruments

12	<b>Wire Tensioner</b>	Graduated Dynamometric Scale	
-	<b>Open Wrench*</b>	10mm Span	
-	<b>Socket Wrench*</b>	10mm Socket	
-	<b>T-Wrench with Adapter*</b>	Accommodates All Half Pin Sizes	

\* Not shown in diagram

