


creed///edge™

Percutaneous Burrs



 GLWmedical

Percutaneous Burrs


Purpose-built for percutaneous foot and ankle procedures, our single-use sterile burrs enable precise osteotomies, bone reshaping, and corrective surgery through incisions as small as 3 – 5mm. Available in cylindrical and conical profiles, each burr is crafted from hardened surgical-grade stainless steel to deliver controlled, high-torque cutting with exceptional wear resistance. Fully compatible with commonly available surgical power units and 4:1 reducers, they integrate directly into your existing MIS workflow.



Shannon Ø2.0,
WL 10mm, TL 55mm
R5-1020-055S



Shannon Ø2.0,
WL 13mm, TL 70mm
R5-1320-070S




Shannon Ø2.2,
WL 22mm, TL 75mm
R5-2222-075S




Shannon Ø3.0,
WL 30mm, TL 100mm
R5-3030-100S



Shannon Ø3.1,
WL 20mm, TL 70mm
R5-2031-070S



Wedge Ø4.3,
WL 13mm, TL 65mm
R5-1343-065S



Ordering Information

Part #	Description
R5-1020-055S	Creed Edge Burr / Shannon Ø2.0, WL 10mm, TL 55mm
R5-1320-070S	Creed Edge Burr / Shannon Ø2.0, WL 13mm, TL 70mm
R5-2222-075S	Creed Edge Burr / Shannon Ø2.2, WL 22mm, TL 75mm
R5-2031-070S	Creed Edge Burr / Shannon Ø3.1, WL 20mm, TL 70mm
R5-3030-100S	Creed Edge Burr / Shannon Ø3.0, WL 30mm, TL 100mm
R5-1343-065S	Creed Edge Burr / Wedge Ø4.3, WL 13mm, TL 65mm

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only. Implants and Instrument kits for treatment of temporary and permanent disabled individuals (patients).

Product availability is subject to regulatory clearance or approval in each applicable market. Please contact your GLW Medical representative if you have questions about the availability of GLW Medical products in your area.

This operative technique guide is intended for healthcare professionals only. It is provided to illustrate a suggested surgical technique for use with a GLW Medical product. It does not constitute medical advice and is not intended to replace the independent clinical judgment of the treating physician. The surgeon is responsible for determining the appropriate technique and whether the product is suitable for a particular patient. Before use, please refer to the applicable product labeling, including the instructions for use and, where applicable, cleaning and sterilization instructions.

GLW Medical is a trade name of GLW, Inc. Creed Edge and Creed are trademarks or registered trademarks of GLW, Inc. All other trademarks are the property of their respective owners. For patent information, please visit www.glwmed.com/patents/.



Legal Manufacturer:

GLW, Inc.
930 Sylvan Ave, Suite 125
Englewood Cliffs, NJ 07632
USA

To order: custsvc@glwmed.com

www.glwmedical.com



Scan for more
product information